


Product Structure – English Translation

Product: One-Piece Open Ostomy Bag (Adult)



Page 1 – Translation Key

Front View = Front Side

Back View = Rear Side

Non-woven ■ / ■■■■ = Non-woven Fabric / Non-woven Ring

Plastic hook face = Plastic Hook Surface

Plastic sheet = Plastic Sheet

Matte sheet (■■■) = Matte Surface Layer


Note: Single-sided skin-coloured non-woven fabric.

Design = Design Review = Review Approval = Approval


Bag Mold = k04 Hydrocolloid Mold = 111 Non-woven Ring = w02

Drawing No. = 50238 Version = V1.0 Effective Date = 2023-07-01


Manufacturer: Zhejiang Ailebao Medical Technology Co., Ltd.




水胶体



无纺布圈



水胶体不干胶



活性炭

一件式开口造口袋 (成人)				图号	50238
				版本	V1.0
袋体模具	k04	设计		生效日期	2023.07.01
水胶体模具	111	审核		浙江爱乐保医疗科技有限公司	
无纺布圈	w02	批准			

Page 2 – Translation Key

Hydrocolloid = ■■■■

Non-woven Ring = ■■■■

Hydrocolloid Adhesive = ■■■■■■

Activated Carbon = ■■■■

Product: One-Piece Open Ostomy Bag (Adult)

Drawing No. 50238 | Version V1.0 | Effective Date 2023-07-01

Manufacturer: Zhejiang Ailebao Medical Technology Co., Ltd.